

Piranha Salt Cell Chlorinator Cleaning Instructions



Step 1 – Remove Salt Cell Chlorinator



***Step 2 – Place Salt Cell Chlorinator
in plastic bag***



***Step 3 – Add water until chlorinator is
fully submerged. Do not exceed more than
1 gallon water per 32 oz. Piranha Salt Cell
Cleaner.***



***Step 4 – Add 1- 32 oz. bottle Piranha
Salt Cell Cleaner to water***



***Step 5 – Soak in solution for 15 – 30
Minutes until clean***



***Step 6 – Remove from bag, flush
contents to drain, dispose of bag and
re-install***